

PPP4R1L Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP07372
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human PPP4R1L. AA range:191-240
Mol wt	49554
Species reactivity	Human, Monkey
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	PPP4R1L Antibody
Synonyms	PPP4R1L; C20orf192; Putative serine/threonine-protein phosphatase 4 regulatory subunit 1-like

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

PPP4R1L (Protein Phosphatase 4 Regulatory Subunit 1 Like (Pseudogene)) is a Pseudogene.


Recommended Dilution

WB: 1: 500 - 1: 2000


ELISA: 1: 10000

Not yet tested in other applications.


Images



Western blot analysis of lysates from COLO, HUVEC, HepG2, Jurkat, and COS cells, using PPP4R1L Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using PPP4R1L antibody.



Western Blot analysis of various cells using PPP4R1L Polyclonal Antibody

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com