

## Rab 24 Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP07536
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from Rab 24 . at AA range: 90-170
Mol wt	23124
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	Rab 24 Antibody
Synonyms	RAB24; Ras-related protein Rab-24

**This product is for research use only, not for use in human, therapeutic or diagnostic procedure.**

### Background

RAB24 (RAB24, member RAS oncogene family) is a small GTPase of the Rab subfamily of Ras-related proteins that regulate intracellular protein trafficking.

### Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 10000

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)