

Rab L5 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP07556
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human RABL5. AA range:128-177
Mol wt	20835
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	Rab L5 Antibody
Synonyms	RABL5; Rab-like protein 5

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

May enhance cellular proliferation. May reduce growth inhibitory activity of CDKN2A.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 5000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com