

Chk2(Phospho Thr432) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP07748
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human Chk2 (Phospho Thr432)
Mol wt	59730
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Chk2
Synonyms	Chk2 ;Phospho Thr432; Serine/threonine-protein kinase Chk2; EC 2.7.11.1; CHK2 checkpoint homolog; Cds1 homolog; Hucds1; hCds1; Checkpoint kinase 2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Catalytic activity:ATP + a protein = ADP + a phosphoprotein.[cofactor:Magnesium].[disease:Defects in CHEK2 are associated with Li-Fraumeni syndrome 2 (LFS2)]

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com