

Ribosomal Protein L11 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP07847
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from Ribosomal Protein L11 . at AA range: 100-180
Mol wt	20252
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Ribosomal Protein L11 Antibody
Synonyms	RPL11; 60S ribosomal protein L11; CLL-associated antigen KW-12

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. RPL11 (ribosomal protein L11) encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L5P family of ribosomal proteins. It is located in the cytoplasm. The protein probably associates with the 5S rRNA. Alternatively spliced transcript variants encoding different isoforms have been found for RPL11. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of RPL11 dispersed through the genome.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 10000

Not yet tested in other applications.

Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4°over night

Storage

-20°C for one year