

CNDP2 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP07896
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Cytosolic non-specific dipeptidase
Synonyms	Cytosolic non-specific dipeptidase ;EC 3.4.13.18;CNDP dipeptidase 2;Glutamate carboxypeptidase-like protein 1;Peptidase A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

CNDP2, also known as tissue carnosinase and peptidase A (EC 3.4.13.18), is a nonspecific dipeptidase rather than a selective carnosinase (Teufel et al., 2003)

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com