

RIOK1 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP07902
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from the Internal region of human RIOK1.
Mol wt	65583
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	RIOK1 Antibody
Synonyms	RIOK1; Serine/threonine-protein kinase RIO1; RIO kinase 1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

RIOK1 (RIO kinase 1) includes two alternatively spliced transcript variants, which encode different isoforms. The function of RIOK1 has not been determined.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 40000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com