

## CL029 Rabbit Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP07984
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human CL029
Mol wt	35750
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Uncharacterized protein C12orf29
Synonyms	Uncharacterized protein C12orf29

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

#### Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for 1 year