

Cleaved-Ephrin-A4 (S170) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08031
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human Cleaved-Ephrin-A4 (S170)
Mol wt	22110
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Cleaved-Ephrin-A4
Synonyms	Cleaved-Ephrin-A4 ;S170; Ephrin-A4; EPH-related receptor tyrosine kinase ligand 4; LERK-4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

May play a role in the interaction between activated B lymphocytes and dendritic cells in tonsils. |Belongs to the ephrin family. |tissue specificity: Expressed in the adult spleen| lymph node| prostate| ovary| small intestine| and colon| and in fetal heart| lung| liver and kidney. Also detected in hematopoietic cell lines. |

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com