

TM4SF19 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08155
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human TM4SF19
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	TM4SF19
Synonyms	TM4SF19 ; Transmembrane 4 L6 family member 19; Osteoclast maturation-associated gene 4 protein; Tetraspan membrane protein OCTM4)

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Belongs to the L6 tetraspanin family. |


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images



Western blot analysis of various lysates, primary antibody was diluted at 1:1000, 4°C overnight, secondary antibody was diluted at 1:10000, 37°C 1hour.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com