

Tmem72 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08193
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human Tmem72 Polyclonal
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Transmembrane protein 72
Synonyms	Transmembrane protein 72 (Kidney-specific secretory protein of 37 kDa); Transmembrane protein 72; Kidney-specific secretory protein of 37 kDa)

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a transmembrane protein which may be expressed specifically in the kidney.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 10000 - 1: 20000

Not yet tested in other applications.

Images



Western blot analysis of various lysate, antibody was diluted at 1000. Secondary antibody was diluted at 1:20000

Storage

-20°C for 1 year