

# Collagen IV alpha 1 (Cleaved-Gly173) Rabbit Polyclonal Antibody

## Description

<b>Product type</b>	Primary Antibody
<b>Code</b>	BT-AP08242
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Size</b>	100ul, 50ul, 20ul
<b>Immunogen</b>	Synthesized peptide derived from human Collagen IV $\alpha$ 1 (Cleaved-Gly173)
<b>Mol wt</b>	183590
<b>Species reactivity</b>	Human, Mouse
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	WB, ELISA
<b>Concentration</b>	1 mg/ml
<b>Full name</b>	Collagen IV alpha 1
<b>Synonyms</b>	Collagen IV $\alpha$ 1 ;Cleaved-Gly173; Collagen alpha-1;IV chain [Cleaved into: Arresten;

**This product is for research use only, not for use in human, therapeutic or diagnostic procedure.**

## Background

disease:Defects in COL4A1 are a cause of brain small vessel disease with hemorrhage

## Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

## Images

No images.

## Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)