

SFT2D3 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08249
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human SFT2D3. AA range:100-150
Mol wt	21790
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	SFT2D3 Antibody
Synonyms	SFT2D3; Vesicle transport protein SFT2C; SFT2 domain-containing protein 3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

SFT2D3 belongs to the SFT2 family. SFT2D3 may be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 40000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com