

Collagen XI alpha 1 (Cleaved-Ala1563) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08255
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human Collagen XI α 1 (Cleaved-Ala1563)
Mol wt	198660
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Collagen XI alpha 1
Synonyms	Collagen XI α 1 ;Cleaved-Ala1563; Collagen alpha-1;XI chain

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

alternative products:Additional isoforms seem to exist. There is alternative usage of exon IIA or exon IIB. Transcripts containing exon IIA or IIB are present in cartilage but exon IIB is preferentially utilized in transcripts from tendon|disease:Defects in COL11A1 are the cause of Marshall syndrome

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com