

COQ5 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08270
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human COQ5
Mol wt	35970
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	2-methoxy-6-polyprenyl-1 4-benzoquinol methylase, mitochondrial
Synonyms	2-methoxy-6-polyprenyl-1; 4-benzoquinol methylase; mitochondrial ;EC 2.1.1.201;Ubiquinone biosynthesis methyltransferase COQ5

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year