

SH2D4A Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08272
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human SH2D4A. AA range:180-220
Mol wt	52727
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	SH2D4A Antibody
Synonyms	SH2D4A; SH2A; SH2 domain-containing protein 4A; Protein SH(2)A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

SH2D4A (SH2 Domain Containing 4A) is a Protein Coding gene. GO annotations related to this gene include phosphatase binding. An important paralog of this gene is HSH2D.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com