

CPNE4 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08311
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human CPNE4
Mol wt	61270
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	CPNE4
Synonyms	CPNE4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


This gene belongs to the highly conserved copine family. It encodes a calcium-dependent| phospholipid-binding protein| which may be involved in membrane trafficking| mitogenesis and development. Alternative splicing results in multiple transcript variants.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com