

SLC10A4 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08350
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from SLC10A4 . at AA range: 340-420
Mol wt	46504
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	SLC10A4 Antibody
Synonyms	SLC10A4; Sodium/bile acid cotransporter 4; Na(+)/bile acid cotransporter 4; Solute carrier family 10 member 4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

SLC10A4 (Solute Carrier Family 10 Member 4) is a Protein Coding gene. GO annotations related to this gene include bile acid: sodium symporter activity. An important paralog of this gene is SLC10A6.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 10000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com