

Sox-17 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08465
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from Sox-17 . at AA range: 90-170
Mol wt	44117
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Sox-17 Antibody
Synonyms	SOX17; Transcription factor SOX-17

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

SOX17 encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The transcription factor SOX-17 may act as a transcriptional regulator after forming a protein complex with other proteins.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 20000

Not yet tested in other applications.

Images



Western Blot analysis of extracts from Jurkat cells, using Sox-17 Polyclonal Antibody. Secondary antibody was diluted at 1:20000 cells nucleus.

Storage

-20°C for one year