

DMRTA Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08530
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 170-250
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Doublesex- and mab-3-related transcription factor A1
Synonyms	Doublesex- and mab-3-related transcription factor A1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

May be involved in sexual development.,Belongs to the DMRT family.,Contains 1 DM DNA-binding domain.,tissue specificity:Expressed in liver, kidney, pancreas, prostate and weakly detected in testis and ovary.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com