

CSMD2 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08549
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human CSMD2
Mol wt	383570
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	IHC-p, IF
Concentration	1 mg/ml
Full name	CSMD2
Synonyms	CSMD2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

IHC-p: 1: 50 - 1: 200

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year