

SUHW4 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08677
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from the Internal region of human SUHW4.
Mol wt	109299
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	SUHW4 Antibody
Synonyms	ZNF280D; KIAA1584; SUHW4; ZNF634; Zinc finger protein 280D; Suppressor of hairy wing homolog 4; Zinc finger protein 634

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

ZNF280D (Zinc Finger Protein 280D) is a Protein Coding gene. An important paralog of this gene is POGZ.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 10000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com