

TAB182 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08782
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human TNKS1BP1. AA range:1601-1650
Mol wt	181782
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	TAB182 Antibody
Synonyms	TNKS1BP1; KIAA1741; TAB182; 182 kDa tankyrase-1-binding protein

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

TNKS1BP1 (Tankyrase 1 Binding Protein 1) is a Protein Coding gene. Among its related pathways are Gene Expression and Regulation of TP53 Activity. Gene Ontology (GO) annotations related to this gene include enzyme binding and ankyrin binding.


Recommended Dilution

WB: 1: 500 - 1: 2000


ELISA: 1: 40000

Not yet tested in other applications.

Images



Western blot analysis of lysates from 293 cells, using TNKS1BP1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using TAB182 Polyclonal Antibody diluted at 1:1000

Storage

-20°C for one year

