

Cyclin B1(Phospho-Ser147) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08828
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phosho peptide around human Cyclin B1 (Ser147)
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Cyclin B1
Synonyms	Cyclin B1 ;Ser147; G2/mitotic-specific cyclin-B1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by this gene is a regulatory protein involved in mitosis. The gene product complexes with p34(cdc2) to form the maturation-promoting factor (MPF). Two alternative transcripts have been found| a constitutively expressed transcript and a cell cycle-regulated transcript| that is expressed predominantly during G2/M phase. The different transcripts result from the use of alternate transcription initiation sites.

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com