

Cyclin D1(Phospho Ser90) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08829
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human Cyclin D1 (Phospho Ser90)
Mol wt	32450
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Cyclin D1
Synonyms	Cyclin D1 ;Phospho Ser90; G1/S-specific cyclin-D1; B-cell lymphoma 1 protein; BCL-1; BCL-1 oncogene; PRAD1 oncogene

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

disease:A chromosomal aberration involving CCND1 may be a cause of B-lymphocytic malignancy| particularly mantle-cell lymphoma (MCL). Translocation t(11;14)(q13;q32) with immunoglobulin gene regions. Activation of CCND1 may be oncogenic by directly altering progression through the cell cycle.|disease:A chromosomal aberration involving CCND1 may be a cause of multiple myeloma

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com