

TCEAL6 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP08884
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from TCEAL6 . at AA range: 1-80
Mol wt	32584
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	TCEAL6 Antibody
Synonyms	TCEAL6; Transcription elongation factor A protein-like 6; TCEA-like protein 6; Transcription elongation factor S-II protein-like 6

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

TCEAL6 (Transcription Elongation Factor A Like 6) is a Protein Coding gene. An important paralog of this gene is TCEAL4.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 10000

Not yet tested in other applications.

Images



Western blot analysis of CoLo using TCEAL6 antibody. Secondary antibody was diluted at 1:20000 cells nucleus.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com