

Thymidine Kinase Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09007
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human KITH. AA range:161-210
Mol wt	25469
Species reactivity	Human, Monkey
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Thymidine Kinase Antibody
Synonyms	TK1; Thymidine kinase; cytosolic

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

TK1 (Thymidine kinase, cytosolic) belongs to the thymidine kinase family and is a regulatory enzyme that can interchange between its dimeric and tetrameric forms, which have low and high catalytic efficiencies, respectively. TK1 is involved in DNA synthesis and an important proliferation marker (PMID: 21178264).


Recommended Dilution

WB: 1: 500 - 1: 2000


ELISA: 1: 5000

Not yet tested in other applications.

Images



Western blot analysis of lysates from HeLa and COS7 cells, using KITH Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from RAW264.7 cells using KITH antibody.

Storage

-20°C for one year

