

DDAH2 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09075
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human DDAH2
Mol wt	31350
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	N
Synonyms	N;G; N;G-dimethylarginine dimethylaminohydrolase 2 ;DDAH-2;Dimethylarginine dimethylaminohydrolase 2;EC 3.5.3.18;DDAHII;Dimethylargininase-2;Protein G6a;S-phase protein

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com