

Desmoyokin(Phospho Ser5782) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09144
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human Desmoyokin (Phospho Ser5782)
Mol wt	647900
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Desmoyokin
Synonyms	Desmoyokin ;Phospho Ser5782; Neuroblast differentiation-associated protein AHNAK; Desmoyokin

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

May be required for neuronal cell differentiation. | Contains 1 PDZ (DHR) domain. | subunit: Interacts with DYSF; the interaction is direct and Ca(2+)-independent. |

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com