

Doublecortin(Phospho-Ser297) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09322
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized phosho peptide around human Doublecortin (Ser297)
Mol wt	N/A
Species reactivity	Human, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Doublecortin
Synonyms	Doublecortin (Ser297); Neuronal migration protein doublecortin; Dublin; Lissencephalin-X; Lis-X)

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


This gene encodes a member of the doublecortin family. The protein encoded by this gene is a cytoplasmic protein and contains two doublecortin domains| which bind microtubules. In the developing cortex| cortical neurons must migrate over long distances to reach the site of their final differentiation. The encoded protein appears to direct neuronal migration by regulating the organization and stability of microtubules. In addition| the encoded protein interacts with LIS1| the regulatory gamma subunit of platelet activating factor acetylhydrolase| and this interaction is important to proper microtubule function in the developing cortex. Mutations in this gene cause abnormal migration of neurons during development and disrupt the layering of the cortex| leading to epilepsy| mental retardation| subcortical band heterotopia ("double cortex" syndrome) in females and lissencephaly ("smooth brain&quo

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year