

DYH2 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09376
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human DYH2
Mol wt	486970
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	IHC-p, IF
Concentration	1 mg/ml
Full name	DYH2
Synonyms	DYH2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Dyneins are microtubule-associated motor protein complexes composed of several heavy| light| and intermediate chains. The axonemal dyneins| found in cilia and flagella| are components of the outer and inner dynein arms attached to the peripheral microtubule doublets. DNAH2 is an axonemal inner arm dynein heavy chain (Chapelin et al.| 1997

Recommended Dilution

IHC-p: 1: 50 - 1: 200

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com