

UBTD2 Polyclonal Antibody

Description

| | |
|--------------------------------|---|
| Product type | Primary Antibody |
| Code | BT-AP09396 |
| Host | Rabbit |
| Isotype | IgG |
| Size | 20ul, 50ul, 100ul |
| Immunogen | The antiserum was produced against synthesized peptide derived from human UBTD2. AA range:181-230 |
| Mol wt | 26190 |
| Species reactivity | Human, Mouse |
| Clonality | Polyclonal |
| Recommended application | WB, ELISA |
| Concentration | 1 mg/ml |
| Full name | UBTD2 Antibody |
| Synonyms | UBTD2; DCUBP; SB72; Ubiquitin domain-containing protein 2; Dendritic cell-derived ubiquitin-like protein; DC-UbP; Ubiquitin-like protein SB72 |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

UBTD2 (Ubiquitin Domain Containing 2) is a Protein Coding gene. An important paralog of this gene is UBTD1.


Recommended Dilution

WB: 1: 500 - 1: 2000


ELISA: 1: 40000

Not yet tested in other applications.

Images



Western blot analysis of UBTD2 Antibody. The lane on the right is blocked with the UBTD2 peptide.



Western blot analysis of the lysates from HUVEC cells using UBTD2 antibody.

Storage

-20°C for one year