

UBTD2 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09396
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human UBTD2. AA range:181-230
Mol wt	26190
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	UBTD2 Antibody
Synonyms	UBTD2; DCUBP; SB72; Ubiquitin domain-containing protein 2; Dendritic cell-derived ubiquitin-like protein; DC-UbP; Ubiquitin-like protein SB72

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

UBTD2 (Ubiquitin Domain Containing 2) is a Protein Coding gene. An important paralog of this gene is UBTD1.


Recommended Dilution

WB: 1: 500 - 1: 2000


ELISA: 1: 40000

Not yet tested in other applications.

Images



Western blot analysis of UBTD2 Antibody. The lane on the right is blocked with the UBTD2 peptide.



Western blot analysis of the lysates from HUVEC cells using UBTD2 antibody.

Storage

-20°C for one year