

EPHA7(Phospho Tyr608) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09663
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human EPHA7 (Phospho Tyr608)
Mol wt	109780
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	EPHA7
Synonyms	EPHA7 ;Phospho Tyr608; Ephrin type-A receptor 7; EC 2.7.10.1; EPH homology kinase 3; EHK-3; EPH-like kinase 11; EK11; hEK11

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Catalytic activity:ATP + a

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com