

ZN-16 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09692
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF292. AA range:1311-1360
Mol wt	218446
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	ZN-16 Antibody
Synonyms	ZNF292; KIAA0530; Zinc finger protein 292

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

May be involved in transcriptional regulation.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 40000

Not yet tested in other applications.

Images



Western blot analysis of ZNF292 Antibody. The lane on the right is blocked with the ZNF292 peptide.

Storage

-20°C for one year