

ZNF575 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09731
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF575. AA range:10-59
Mol wt	26763
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	ZNF575 Antibody
Synonyms	ZNF575; Zinc finger protein 575

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

May be involved in transcriptional regulation.


Recommended Dilution

WB: 1: 500 - 1: 2000


ELISA: 1: 5000

Not yet tested in other applications.

Images



Western blot analysis of the lysates from HeLa cells using ZNF575 antibody.



Western blot analysis of ZNF575 Antibody. The lane on the right is blocked with the ZNF575 peptide.

Storage

-20°C for one year