

ZnT-2 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09749
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human SLC30A2. AA range:150-200
Mol wt	35178
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	ZnT-2 Antibody
Synonyms	SLC30A2; ZNT2; Zinc transporter 2; ZnT-2; Solute carrier family 30 member 2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by SLC30A2 (solute carrier family 30 member 2) is a zinc transporter that acts as a homodimer. The encoded protein plays a role in secreting zinc into breast milk. Two transcript variants encoding different isoforms have been found for SLC30A2.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 5000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com