

HEAT1 Polyclonal Antibody

Description

| | |
|--------------------------------|---|
| Product type | Primary Antibody |
| Code | BT-AP09856 |
| Host | Rabbit |
| Isotype | IgG |
| Size | 100ul, 50ul, 20ul |
| Immunogen | Synthesized peptide derived from part region of human protein |
| Mol wt | N/A |
| Species reactivity | Human, Rat, Mouse |
| Clonality | Polyclonal |
| Recommended application | WB, ELISA |
| Concentration | 1 mg/ml |
| Full name | HEAT repeat-containing protein 1 |
| Synonyms | HEAT repeat-containing protein 1 ;Protein BAP28 |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Involved in nucleolar processing of pre-18S ribosomal RNA. Involved in ribosome biosynthesis.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,Belongs to the HEATR1/UTP10 family.,Contains 1 HEAT repeat.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com