

FA11 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09939
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human FA11 (heavy chain Cleaved-Arg387)
Mol wt	68750
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	FA11
Synonyms	FA11 ;heavy chain; Cleaved-Arg387; Coagulation factor XI; FXI; EC 3.4.21.27; Plasma thromboplastin antecedent; PTA; Coagulation factor XIa heavy chain; Coagulation factor XIa light chain;

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes coagulation factor XI of the blood coagulation cascade. This protein is present in plasma as a zymogen| which is a unique plasma coagulation enzyme because it exists as a homodimer consisting of two identical polypeptide chains linked by disulfide bonds. During activation of the plasma factor XI| an internal peptide bond is cleaved by factor XIIa (or XII) in each of the two chains| resulting in activated factor XIa| a serine protease composed of two heavy and two light chains held together by disulfide bonds. This activated plasma factor XI triggers the middle phase of the intrinsic pathway of blood coagulation by activating factor IX. Defects in this factor lead to Rosenthal syndrome| a blood coagulation abnormality.

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year