

FA46D Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09949
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human FA46D
Mol wt	42790
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	FA46D
Synonyms	FA46D

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


Antibodies against the protein encoded by this gene were found only in plasma from cancer patients. While it may be a target for immunotherapy| the function of this gene is unknown.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from MDA-MB cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com