

FA7 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09966
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human FA7 (light chain Cleaved-Ala61)
Mol wt	51260
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	FA7
Synonyms	FA7 ;light chain; Cleaved-Ala61; Coagulation factor VII; EC 3.4.21.21; Proconvertin; Serum prothrombin conversion accelerator; SPCA; Eptacog alfa; Factor VII light chain; Factor VII heavy chain;

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Catalytic activity: Selective cleavage of Arg-Ile bond in factor X to form factor Xa. | disease: Defects in F7 are the cause of factor VII deficiency

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com