

FEZ2 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10025
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human FEZ2
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	FEZ2
Synonyms	FEZ2; Fasciculation and elongation protein zeta-2; Zygin II; Zygin-2)

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


This gene is an ortholog of the *C. elegans unc-76* gene| which is necessary for normal axonal bundling and elongation within axon bundles. Other orthologs include the rat gene that encodes zygin II| which can bind to synaptotagmin.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com