

FGF19 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10038
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human FGF19
Mol wt	23760
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	FGF19
Synonyms	FGF19; Fibroblast growth factor 19; FGF-19

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

May be involved in brain development during embryogenesis. [miscellaneous:Contrarily to other members of the family that can bind several FGF receptors FGF19 is specific for FGFR4. |Belongs to the heparin-binding growth factors family. |tissue specificity:Expressed in fetal brain| cartilage| retina| and adult gall bladder.]

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com