

FoxO3A (Acetyl Lys569) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10072
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human FoxO3A (Acetyl Lys569)
Mol wt	74030
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	FoxO3A
Synonyms	FoxO3A ;Acetyl Lys569; Forkhead box protein O3; AF6q21 protein; Forkhead in rhabdomyosarcoma-like 1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

disease:A chromosomal aberration involving FOXO3 is found in secondary acute leukemias. Translocation t(6;11)(q21;q23) with MLL/HRX.[Transcriptional activator which triggers apoptosis in the absence of survival factors] including neuronal cell death upon oxidative stress. Recognizes and binds to the DNA sequence 5'-

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com