

GAPDH Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10120
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Purified recombinant Full length GAPDH (Yeast)
Mol wt	N/A
Species reactivity	Yeast
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	GAPDH
Synonyms	GAPDH; G3PD; G3PDH; cb609; HEL-S-162eP; MGC102546; MGC103190; MGC103191

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A


Recommended Dilution

WB: 1: 5000 - 1: 20000

ELISA: 1: 20000

Not yet tested in other applications.

Images



Western Blot analysis of Yeast cells using GAPDH rabbit pAb(Yeast). Secondary antibody was diluted at 1:20000

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com