

HS12A Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10162
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 230-310
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Heat shock 70 kDa protein 12A
Synonyms	Heat shock 70 kDa protein 12A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Belongs to the heat shock protein 70 family.,tissue specificity:Widely expressed with highest levels in brain, kidney and muscle.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com