

GLT13 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10189
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human GLT13
Mol wt	61160
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Polypeptide N-acetylgalactosaminyltransferase 13
Synonyms	Polypeptide N-acetylgalactosaminyltransferase 13 ;EC 2.4.1.41;Polypeptide GalNAc transferase 13;GalNAc-T13;pp-GaNTase 13;Protein-UDP acetylgalactosaminyltransferase 13;UDP-GalNAc:polypeptide

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com