

## Granulins Rabbit Polyclonal Antibody

### Description

<b>Product type</b>	Primary Antibody
<b>Code</b>	BT-AP10278
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Size</b>	100ul, 50ul, 20ul
<b>Immunogen</b>	Synthesized peptide derived from human Granulins
<b>Mol wt</b>	65230
<b>Species reactivity</b>	Human, Rat, Mouse
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	WB, ELISA
<b>Concentration</b>	1 mg/ml
<b>Full name</b>	Granulins
<b>Synonyms</b>	Granulins; Granulins; Proepithelin; PEPI; Acrogranin; Paragranulin; Granulin-1; Granulin G; Granulin-2; Granulin F; Granulin-3; Granulin B; Granulin-4; Granulin A; Granulin-5; Granulin C; Granulin-6; Granulin D; Granulin-7; Granulin E;

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

disease:Defects in GRN are the cause of ubiquitin-positive frontotemporal dementia (UP-FTD)

### Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for 1 year