

GTPBA Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10316
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human GTPBA
Mol wt	42570
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	GTPBA
Synonyms	GTPBA

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


Small G proteins| such as GTPBP10| act as molecular switches that play crucial roles in the regulation of fundamental cellular processes such as protein synthesis| nuclear transport| membrane trafficking| and signal transduction (Hirano et al.| 2006

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com