

H2B2C Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10338
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human H2B2C
Mol wt	21230
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Putative histone H2B type 2-C
Synonyms	Putative histone H2B type 2-C ;Histone H2B.t;H2B/t

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year